

CHART PROJECTOR

CP-500

CHART PROJECTOR CP-500



SHIN-NIPPON by Rexxam

by Rexxam CHART PROJECTOR

50charts and various masks

Landolt, Letter, Child, Snellen, 50 kinds of charts, and Horizontal line, Vertical line, one chart, various masks can be selected by the remote control box.

Fast, slient chart disk rotation

It moves to the following chart in a maximum of 0.3 seconds.

User Program mode

2 kinds of programs can be selected. 48 kinds (maximum) of chart and mask are memorizable.

Room light swicth


By using optional deviator, operation of a room light is available by the remote control.


- 1 Lamp switch 4 Chart/Mask switch
- 2 Room light switch
- 5 Horizontal line switch 8 Back switch
- One chart switch 10 Eyeball chart switch

- 3 Program switch (1,2) 6 Vertical line switch
 - 9 Next switch
- **SPECIFICATIONS** Number of chrats and masks 50charts and 38masks min.80-max.320Cd/m2(200Cd/m2 at 5m) Luminance Chart and mask selection speed 0.3seconds 2 kinds(48charts or masks) Program ±10° Tilting angle Auto power off After 5 min. 12V, 30W halogen lamp Lamp **Projection distance** 2~7m Projection magnification 30 x (at 5m) Projection size 335mm x 222mm (at 5m) Dimentions W248xL336xH234mm Weight 4.5kg 100~240V Power source Power consumption 150VA STANDARD ACCESSORIES

REMOTE CONTROL BOX, SUM4 BATTERY X 2, SCREEN, TEST SCALE SHEET, DUST COVER, SPARE LAMP, SPARE FUSE **OPTIONAL ACCESSORIES**


TABLE BASE STAND, WALL MOUNT, CHART PROJECTOR ARM, POLA FILTER GLASS, R&G GLASS, DEVIATOR







Contact MEC Sales Division 2-8-4, Kandatsukasa-machi, Chiyoda-ku Tokyo, 101-0048, Japan TEL 81-3-3256-7701 FAX 81-3-3256-7702 E-mail: eye@rexxam.co.jp URL: http://www.rexxam.co.jp URL: http://www.shin-nippon.jp



Adjustment of the brightness

The brightness of the projector can be changeable by the control switch according to the brightness of the operation room.



Manufacture



TOWA MEDICAL INSTRUMENTS Co., Ltd. 1698-2 Oaza Nakano Nakano-shi, Nagano-ken, 383-0013 Japan